

## PDC-VGC VACUUM GAUGE & DIGITAL METER



PDC-VCG and accessories



PDC-VCG shown with High Power Expanded Plasma Cleaner (sold separately)

## PDC-VCG (115V) PDC-VCG-2 (230V)

## **FEATURES**

- Thermocouple vacuum gauge sensor
- Digital meter displays pressure range of 1 mTorr to 760 Torr
- Tee adapter included to connect vacuum gauge to plasma cleaner front door
- 1/4" OD hose adapter included to connect input gas line to metering valve on plasma cleaner
- Weight: 5 lbs | Size: 5" H x 10" W x 8" D

## COMPATIBLE MODELS

- PDC-32G / PDC-32G 2 Basic Plasma Cleaner
- PDC-001 / PDC-002 Expanded Plasma Cleaner
- PDC-001-HP / PDC-002-HP High Power Expanded Plasma Cleaner